

Medlink's Mini PACS (Picture Archive and Communication System) allows for a cost effective transition from film to digital for Private Practices such as podiatrists, chiropractors and medical practices. Using the Orex ACL CR systems and Cedara's I-View software gives the practice a cost effective Mini PACS. Medlink's platform is designed to utilize commercial off the shelf hardware and operating systems for easy repair and user friendly navigation. Practitioners no longer have to pay for film storage, darkrooms, water and chemistry. The digital conversion allows practices to be more competitive in an environment of increased cost and reduced reimbursements.

Typical Configuration:

- Orex ACL2 Computed Radiography System
- Dell Dimension 4700 PC
- Cedara HIS / RIS HL 7 Patient, Study Update / Merge (IHE Compliant) (Option)
- Dell PowerEdge 2800 Server - 300 GB Hard drive - stores 300,000 images
- Cedara I-View Software v6.x or Cedara I-Reach "Single Practitioner" Web Server/Client/Archive Module for A Single User



## Typical Private Practic Configuration

