

Mobile Hybrid CT
DeteCT series

JPI Healthcare Solutions

3555 Veterans Memorial Hwy, Unit D, Ronkonkoma NY 11779

sales@jpihealthcare.com



www.jpianimalhealth.com/detect-vet

Copyright (c) JPI Healthcare Solutions, Inc. All Rights Reserved.



X-ray Grids



Imaging System



Detector




Software

Innovative Leader
In X-ray Imaging Solution

VETERINARY
Imaging Solution

Mobile Hybrid CT
DeteCT Vet



Hybrid (3-in-1)



Mobility



Usability






Compact (Slim)

www.jpianimalhealth.com


High-definition CT system that provides accurate 3D images with AI-powered acquisition software

Mobile Hybrid CT

DeteCT Vet




Dental Imaging



Spine Imaging


CT



TOMOGRAPHY

- Fast scan time (4-20 sec)
- MPR viewer
- Superior image quality
- Fast reconstruction time (25 sec)


FLO



FLUOROSCOPY

- Large FoV: 17" x 17"
- Frame rate: 10 ~ 30 fps
- Real-time image processing
- ABC (Auto Brightness control)
- DSA (Digital Subtraction Angiography)


RAD



RADIOGRAPHY

- Simple workflow
- Image preview
- Advanced image processing
- Powerful measurement tools
- Fast image processing time (< 0.8 sec)

Flexibility and Efficiency do not have to come at the cost of performance compromise.



Control Touch Panel



User-friendly and Intuitive UX and UI paired with easy maneuverability.

Customer Service

Our outstanding customer support team is consisted of engineers with experience and knowledge. Also, we supply remote training session for our valued partners so they may enhance their understanding and supply the best service possible to end users.

