



# High Frequency Portable X-Ray 9020H.F



# 9020

High Frequency Portable Generator

## Specifications

---

- Ultra light weight-13lbs
  - 40-90kV step /20mA, 18mA, 16mA
  - 1.60 kW Max Output
  - 100 kHz h.f switching inverter
  - 2 Year Warranty
  
  - Dual laser pointers (26") & measuring tape
  - Built in collimator
  - Reversible digital display of kV & mAs
  - Soft touch user panel
  - Overload protection for kV,mA time and temp.
  - Electronic voltage compensation
  - Self manageable APR functions with 5 memory keys
  - Foam grips on each side of handle for easy handling
  - Laser pointers and collimator light activate with a touch of the exposure switch
  - Carrying case for easy transport
  - Ideal for equine applications
- 

Authorized Dealer

Raytech Diagnostics

613-799-4171

<https://raytechdiagnostics.com/>