



TruePortable

Versatile Imaging. On-the-go.

Designed with veterinary practitioners in mind, TruePortable™ stands out as a superior option—lighter, more maneuverable, and offering enhanced protection compared to its competitors. Easily navigable on stairs, in tight spaces, and across various field conditions, TruePortable™ is tailored to meet the unique needs of veterinary professionals.



Our Solutions

With superior image quality and intuitive veterinary imaging software, the ExamVue™ Duo digital imaging system delivers fast image acquisition, easy workflow management, and lower dose for your patients.

The ExamVue™ Duo system comes equipped with a strong suite of veterinary tools, and is available in any configuration for your practice needs.

ExamVue™ Duo Package



ExamVue Duo
14"x17" Wireless



Laptop
Workstation



X-ray Grid



ExamVue Duo Veterinary
Acquisition Software



DICOM
Modality Worklist

Choose Your PACS



LAN PACS



Cloud PACS

Powerful & Versatile

- ▶ 82lbs/37kg (Generator and Stand)
- ▶ 140µm Csl panel
- ▶ 100 kV generator
- ▶ Control the generator directly from the laptop

Unmatched Coverage Program

- ▶ Up to 5 years of PANEL WARRANTY
- ▶ Up to 5 years of REMOTE SOFTWARE SUPPORT
- ▶ Up to 5 years of DROP COVERAGE PROGRAM
- ▶ 3 years of LAPTOP WORKSTATION WARRANTY
- ▶ 2 years of X-RAY GENERATOR WARRANTY
- ▶ 1 year of MOBILE STAND WARRANTY






Premium Protection

TruePortable™ combines portability, durability, power, and advanced image processing to provide exceptional images. We built TruePortable™ with protection of your investment in mind. Protective encasements and carry bag protect your detector and laptop during transportation and use.

Great protection will allow your team to go about their day knowing that your most valuable components are well-protected no matter where you take them.



ExamVue™ Duo Veterinary Acquisition Software

Utilizing cutting-edge and user-friendly imaging tools, our ExamVue™ Duo veterinary acquisition software effortlessly generates clear, concise, and high-quality images, proving effective even with older generators. Offering reliability and convenience, the ExamVue™ Duo seamlessly integrates into any veterinary practice.



ExamVue™ PACS includes built-in veterinary diagnostic tools and can be enhanced with the JPI Cloud PACS. This cloud-based solution extends your image accessibility to any internet-connected location, facilitating seamless sharing with veterinary professionals and patients within your practice or globally.

Ergonomic Design



Folded Vertically

Folded Horizontally

Unfold and be ready to take an X-ray immediately

User-Friendly Design

Simple & Fast Workflow

Convenient Diagnosis Function

Contact us:

3555 Veterans Memorial Hwy, Unit D
Ronkonkoma, NY 11779
sales@jpihealthcare.com
(516) 513-1330 Ext. 3



www.jpihealthcare.com