



ExamvueCR Pro Readers



Key Benefits

- Affordable - Fast Speed
- Wall-mountable
- Compact - Portable
- Crystal-Clear Image Quality
- Easy-to-use software
- Light weight at 42 lbs
- 3 year parts warranty



Mobile application

Technical Specifications

- 3 Cassette sizes: 14x 17, 10 x 12 and 8 x 10 in
- 70 large plates/hour
- Weight: 42 lbs
- Size: 5.4 x 18 x 31.6 in
- Auto-alignment after transportation



Veterinary Dental Cassette



Key Benefits

- Cost-Effective
- Scan Dental IPs without Dedicated Dental Scanner
 - Note resolution 100 microns



Technical Specifications

- Cassette size 10x 12IP sizes: Package comes with two #2 and one #4
- Sizes available :0, 1, 2, 3, 4c and 4




Acquisition



Key Benefits

- Simple and intuitive
- Customizable system
- Embedded viewer for measurements and annotation
- Dicom Query/Send/Receive, Dicom print enabled
- CD/DVD Burn with embedded viewer included

Intuitive Workflow



1 2
3 4



RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com
613-799-4171