



Specifications	Value	Specifications	Value
<b>X-ray Generator</b>		<b>Mobile Stand</b>	
CubeX 28 HF X-ray Generator	2.7 kW, 40 - 100 kVp, 40 mA Max 28 lbs (13 kg)	Weight	92 lbs (41.7 kg) with generator and cables
Controls	Integrated controls with ExamVue Duo Software Inverted Control panel for use on a stand or table	Maximum SID to Floor	79" (200 cm)
Panel		Large Wheels	10" (25.4 cm)
ExamVue Duo EV14W	14" x 17" Wireless Detector Assembled in the US 140 $\mu$ m pixel pitch, 3072 x 2560 pixels active area  Direct Deposit CSI with AED  Dual Band Wireless Support (2.4 Ghz and 5 Ghz)  3.4 kg (7.5 lbs)  Quick connect magnetic cable Includes 2 batteries and battery charger	Folded Dimensions	27" (61 cm) x 57" (137 cm) x 27" (61 cm)
Panel Protection Options		Integrated Cabling	Exposure Switch Extension
Mobile DR Panel Holder	14" x 17" with Single Handle Long Edge	<b>Acquisition Software</b>	
Full Panel Encasement With Snap-On Grid	103 LPI 14" x 17" with Aluminum Frame for Mobile	ExamVue Duo	Acquisition Software for Medicine
		Acquisition Laptop	* Specs subject to change without notice  Windows 10 PRO  16 GB RAM, 512 GB SSD Hard Drive  i5 10th GEN Processor  15.6" display
		<b>Support</b>	
		Remote Service Support Software	Team Viewer Remote Access Software Included (allows for WEB-Based Remote dial-up support)

## Contact us:

RAYTECH-DIAGNOSTICS  
613-799-4171  
sales@raytechdiagnostics.com  
<https://raytechdiagnostics.com>