

ScanX® Duo

Acquire | Store | View | Share

The Most Complete Dental CR Available



Unmatched Workflow Advantages

The ScanX Duo Digital Imaging System with In-Line Erase is another innovative imaging product from ALLPRO Imaging, a leading manufacturer of veterinary equipment since 1962. The ScanX Duo is designed to process intraoral imaging plates. The device has been designed and manufactured using state-of-the-art technology to give many years of dependable service.

For the busy veterinary practice, the ScanX Duo offers workflow advantages by being able to scan all plate sizes, 0, 1, 2, 3, 4, 5 and 6. In addition, with multiple processing slots, two plates can be scanned simultaneously. And since the patient hygiene is always a concern, the ScanX Duo scans without cross contamination from barrier to device.

To complement your ScanX Duo, you need imaging software to acquire, store, view and share images. TigerView, a Televere Systems product, will significantly improve your work flow for better patient management and diagnosis. Quickly and easily access all digital image files and share them with colleagues and patients. Window leveling and viewer enhancement tools allow for complete control of image display.



21 LBS

21 lbs.



Delivers 22 line pairs/mm and 1,100 dots per inch



Enables staff to complete revenue generating procedures for your practice



Detect pathologies under the gum




Capable of reading images simultaneously from 2 PSP Imaging Plates



Size 6 and 4 plates allow for shooting films in the OR on extremities

ScanX® Duo





The Most Complete Dental CR Available



SPECIFICATIONS

Part Number	PN D1000V
Dimensions	10.375" W x 10.0" L x 12.0" H (26.352 cm x 25.4cm x 30.48 cm)
Resolution	Up to 22 lp/mm
Weight	21 lbs (9.5 kg)
Electrical	100 to 240VAC +/- 10%, 50/60Hz
Laser Classification	Class 1 Laser Product Compliance with FDA 21 CFR 1040.10 and IEC 60825-1
Port	USB
Imaging Media	Two processing slots accept Phosphor Storage Plates intraoral sizes 0 through 6
System Includes	Power Supply and Computer Connector Cord
Optional Accessories	PSPs, Intraoral: sizes 0-6 PSP Cleaning Wipes: package of 50 Barrier Envelopes: sizes 0-6 Plate Guides: sizes 0-4 Plate Transfer Box: transports exposed intraoral PSPs ScanX Cleaning Sheets: packages of 12 and 50

- Digital, Wireless and Versatile
- The ScanX PSP Advantage
- Unmatched Detail and Image Quality



[https://www.raytechdiagnostics.com/
sales@raytechdiagnostics.com](https://www.raytechdiagnostics.com/sales@raytechdiagnostics.com)
613-799-4171