

STORED ENERGY GENERATOR – SE 30

- NO BATTERIES!
- 30KW FROM A 240 VAC, 15 AMP RECEPTACLE

Utilizing the latest in high density electrolytic capacitor storage techniques, the SE-30 has no batteries to slowly lose their charge and require expensive replacement. Rated at more than 10,000 hours of continuous usage, these space age capacitors are at the heart of the stored energy solution provided by the SE-30.



Stored Energy High Frequency X-ray Generator offers a “no compromise” solution for customers who have a power problem. The 30 kW SE-30 will provide all the x-ray energy needed for your Veterinary, Chiropractic or Medical Clinics. Plug it in to a 240 VAC, 15 amp outlet and you’re ready to calibrate and take x-rays.



Technical Information:

Output power – 30kW
Technique factors – 400 mA, 125 kVp, 30 mAs (TV) 100 mAs (TH)
mA Range – 25 to 400 mA
kVp Range – 40 to 125 kVp
Time – 2ms to 2 sec.
Input power – 240 VAC, Single Phase, 15 Amp Breaker
Size – 44" x 24" x 12.5"
Weight – 450 lbs (205 kgs)

Features:

APR – optional
Bucky Control – optional
Tube Protection – Yes
Note: If the exact mAs is not selectable, choose the next closest.



RAYTECH DIAGNOSTICS
<https://www.raytechdiagnostics.com>
sales@raytechdiagnostics.com
613-721-1982