

Weight Bearing Protection



Our new Weight Bearing Protective Cover has been designed to accommodate Digital DR Panels and CR Cassettes.

Our covers are strong yet lightweight, offering a high degree of protection during weight bearing procedures.

Our latest cover is constructed of a high density polyethylene frame and a carbon fiber top, offering protection up to 1300 lbs. Our easy lift handle makes removal and carrying a breeze.

The Weight Bearing Cover has been tested by an accredited independent laboratory.



RC Imaging
your radiology resource

Weight Bearing Protection



Smoked Polycarbonate



Carbon Fiber

- Available in Polycarbonate or Carbon Fiber
- Rated up to 1300 lbs.
- Easy lift handle
- Non-skid rubber feet
- Multiple sizes available
- Printed target area

This product is perfect for bariatric, equine, and podiatry procedures.



Easy Lift Handle

Protect your **investment** use a **cover!**



RC Imaging
your radiology resource

RAYTECH DIAGNOSTICS

<https://www.raytechdiagnostics.com/>

sales@raytechdiagnostics.com

613-799-4171

