


Grid Encasements



Avoid future grid expense and encase your grid!
Extend the life of your Grid by using one of our Grid Encasements built to your specifications. They are available in multiple sizes to meet your needs.



RC Imaging
your radiology resource

GridPro Grid Encasements

Polycarbonate with Channels GEC



Extreme Durability • Cassette slides easily in and out of channels
Minimal patient disruption during multiple exposures

Aluminum with Channels GECA



Vinyl coated target area for easy alignment and patient comfort
Reinforced welded corners • Minimal patient disruption during multiple exposures

Aluminum Drop-On Grid Cap GEDA



Vinyl coated target area for easy alignment and patient comfort
Open corners allow easy removal • Cost effective



RC Imaging
your radiology resource

RAYTECH DIAGNOSTICS

<https://www.raytechdiagnostics.com/>
sales@raytechdiagnostics.com 613-799-4171

